

Product Datasheet

Alpha-mannosidase II Antibody (PSH02-67) NBP3-32565

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-32565

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-32565



NBP3-32565**Alpha-mannosidase II Antibody (PSH02-67)**

Product Information	
Unit Size	100 ul
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	PSH02-67
Preservative	0.05% Sodium Azide
Isotype	IgG
Purity	Protein A purified
Buffer	1*TBS (pH7.4), 0.1% BSA and 40% Glycerol
Target Molecular Weight	131 kDa

Product Description	
Description	Novus Biologicals Rabbit Alpha-mannosidase II Antibody (PSH02-67) (NBP3-32565) is a recombinant monoclonal antibody validated for use in IHC, WB and Flow. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	4124
Gene Symbol	MAN2A1
Species	Human, Mouse, Rat
Immunogen	Recombinant protein within human Alpha-mannosidase II aa 51-500 / 1,144. (Uniprot: Q16706)

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Flow Cytometry, Immunohistochemistry
Recommended Dilutions	Western Blot 1:1000, Flow Cytometry 1:1000, Immunohistochemistry, Immunohistochemistry-Paraffin 1:1000



Images

Western Blot: Alpha-mannosidase II Antibody (PSH02-67) [NBP3-32565]
- Western blot analysis of Alpha-mannosidase II on different lysates with Rabbit anti-Alpha-mannosidase II antibody (NBP3-32565) at 1/1,000 dilution.

Lane 1: HepG2 cell lysate
Lane 2: Huh7 cell lysate
Lane 3: JAR cell lysate
Lane 4: U-87 MG cell lysate
Lane 5: HUVEC cell lysate
Lane 6: A375 cell lysate
Lane 7: MDA-MB-468 cell lysate
Lane 8: NCI-H226 cell lysate
Lane 9: Human liver tissue lysate
Lane 10: human lung tissue lysate
Lane 11: human kidney tissue lysate

Lysates/proteins at 30 ug/Lane.

Predicted band size: 131 kDa
Observed band size: 131-151 kDa

Exposure time: 5 minutes;

4-20% SDS-PAGE gel.




Proteins were transferred to a PVDF membrane and blocked with 5% NFDM/TBST for 1 hour at room temperature. The primary antibody (NBP3-32565) at 1/1,000 dilution was used in 5% NFDM/TBST at room temperature for 2 hours. Goat Anti-Rabbit IgG - HRP Secondary Antibody at 1/50,000 dilution was used for 1 hour at room temperature.

Immunohistochemistry: Alpha-mannosidase II Antibody (PSH02-67) [NBP3-32565] - Immunohistochemical analysis of paraffin-embedded human colon cancer tissue with Rabbit anti-Alpha-mannosidase II antibody (NBP3-32565) at 1/1,000 dilution.

The section was pre-treated using heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0) for 20 minutes. The tissues were blocked in 1% BSA for 20 minutes at room temperature, washed with ddH₂O and PBS, and then probed with the primary antibody (NBP3-32565) at 1/1,000 dilution for 1 hour at room temperature. The detection was performed using an HRP conjugated compact polymer system. DAB was used as the chromogen. Tissues were counterstained with hematoxylin and mounted with DPX.

Flow Cytometry: Alpha-mannosidase II Antibody (PSH02-67) [NBP3-32565] - Flow cytometric analysis of U-87 MG cells labeling Alpha-mannosidase II.

Cells were washed twice with cold PBS and resuspend. Then stained with the primary antibody (NBP3-32565, 1/1,000) (red) compared with Mouse IgG Isotype Control (green). After incubation of the primary antibody at +4°C for an hour, the cells were stained with a iFluor™ 488 conjugate-Goat anti-Mouse IgG Secondary antibody at 1/1,000 dilution for 30 minutes at +4°C. Unlabelled sample was used as a control (cells without incubation with primary antibody; black).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950

Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400

Toll Free: 855.668.8722

Fax: 905.827.6402

canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15

Fax: (44) (0) 1235 533420

info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com

Technical Support: nb-technical@bio-techne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

Products Related to NBP3-32565

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-32565

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

