

Product Datasheet

KRAS Antibody (D113) - Mutant - BSA Free NBP3-43582

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-43582

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-43582



NBP3-43582

KRAS Antibody (D113) - Mutant - BSA Free


Product Information	
Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	D113
Preservative	0.02% Proclin 300
Isotype	IgG1
Purity	Protein A purified
Buffer	PBS

Product Description	
Description	Novus Biologicals Mouse KRAS Antibody (D113) - Mutant - BSA Free (NBP3-43526) is a recombinant monoclonal antibody validated for use in ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	3845
Gene Symbol	KRAS
Species	Human
Specificity/Sensitivity	This antibody binds the G12D mutant form of human KRAS protein. This ras protein binds GDP/GTP and possess intrinsic GTPase activity, that results in conversion of nucleotide guanosine triphosphate (GTP) into guanosine diphosphate (GDP). This protein plays an important role in many signal transduction pathways and regulation of cell proliferation. Somatic KRAS mutations are highly observed in leukemias, colorectal cancer, pancreatic cancer and lung cancer. KRAS mutation is G12D which is estimated to be present in up to 37% pancreatic cancers and over 12% of colorectal cancers.
Immunogen	The original antibody was generated by immunizing mice with a synthetic dodecapeptide 'KLVVVGADGVGK'. (Uniprot# P01116-2)

Product Application Details	
Applications	ELISA
Recommended Dilutions	ELISA

Images

Immunocytochemistry/ Immunofluorescence: KRAS Antibody (D113) - Mutant - BSA Free [NBP3-43582] - Immunofluorescence analysis of paraformaldehyde fixed THP-1 cells on coverslips stained with the chimeric version of this antibody, NBP3-43583 (1:100 dilution) for 1h followed by Alexa Fluor® 488 secondary antibody (1:1000 dilution), showing membrane and cytoplasmic staining. The nuclear stain is DAPI (blue). Panels show, from left-right, top-bottom, NBP3-43583, DAPI, merged channels and an isotype control. The isotype control was an unknown specificity antibody followed by staining with Alexa Fluor® 488 secondary antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950

Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400

Toll Free: 855.668.8722

Fax: 905.827.6402

canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15

Fax: (44) (0) 1235 533420

info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com

Technical Support: nb-technical@bio-techne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

Products Related to NBP3-43582

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
H00003845-P01-10ug	Recombinant Human KRAS GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-43582

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

