

# Product Datasheet

## PEX5 Antibody (OT16C4) [DyLight 550] NBP3-44634R

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-44634R](http://www.novusbio.com/NBP3-44634R)

Updated 5/20/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-44634R](http://www.novusbio.com/reviews/destination/NBP3-44634R)



**NBP3-44634R**

PEX5 Antibody (OT16C4) [DyLight 550]

| Product Information |   |
|---------------------|---|
| Unit Size           | 0.1 ml  |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.  |
| Clonality           | Monoclonal  |
| Clone               | OT16C4  |
| Preservative        | 0.05% Sodium Azide  |
| Isotype             | IgG2b   |
| Conjugate           | DyLight 550   |
| Purity              | >95%  |
| Buffer              | 50mM Sodium Borate  |

| Product Description |   |
|---------------------|---|
| Host                | Mouse   |
| Gene ID             | 5830  |
| Gene Symbol         | PEX5  |
| Species             | Human, Mouse, Rat   |
| Immunogen           | Full length human recombinant protein of human PEX5 (NP_000310) produced in HEK293T cell. |
| Notes               | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.         |

| Product Application Details |  |
|-----------------------------|--|
| Applications                | Western Blot, Immunocytochemistry/ Immunofluorescence                  |
| Recommended Dilutions       | Western Blot, Immunocytochemistry/ Immunofluorescence                  |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-44634R**

---

|               |  |
|---------------|--|
| NBP2-27231R   | Mouse IgG2b Isotype Control (MPC-11) [DyLight 550] |
| NBP1-87185PEP | PEX5 Recombinant Protein Antigen                   |
| 8090-ZN-025   | PEX5   |
| NBL1-14304    | PEX5 Overexpression Lysate                         |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-44634R](http://www.novusbio.com/reviews/submit/NBP3-44634R)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

