

Product Datasheet

Bone marrow stromal cell antigen 1/CD157 Antibody (DMC444) - Chimeric NBP3-48341-100ug

Unit Size: 100 ug

Store at -20C in powder form. Store at -80C once reconstituted.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-48341

Updated 8/13/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-48341



NBP3-48341-100ug

Bone marrow stromal cell antigen 1/CD157 Antibody (DMC444) - Chimeric

| Product Information | |
|------------------------------------|---|
| Unit Size | 100 ug |
| Concentration | LYOPH mg/ml |
| Storage | Store at -20C in powder form. Store at -80C once reconstituted. |
| Clonality | Monoclonal |
| Clone | DMC444 |
| Preservative | Trehalose |
| Reconstitution Instructions | Dissolve the lyophilized antibody in sterile and low endotoxin water. The recommend volume of water is lot specific. Please refer to the instructions printed on the product vial or contact technical support. |
| Isotype | IgG1 |
| Purity | Antigen Affinity-purified from cell culture supernatant |
| Buffer | PBS, pH 7.4, 5-8% Trehalose |

| Product Description | |
|----------------------------|---|
| Description | This product is shipped at ambient temperature. Store at -20C to -80C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80C (Avoid repeated freezing and thawing). |
| Host | Human |
| Gene ID | 683 |
| Gene Symbol | BST1 |
| Species | Human |
| Immunogen | Bone marrow stromal cell antigen 1/CD157 recombinant protein expressed in mammalian cells (Accession# Q10588). The exact amino acid sequence is proprietary. |

| Product Application Details | |
|------------------------------------|----------------------|
| Applications | Flow Cytometry |
| Recommended Dilutions | Flow Cytometry 1:100 |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-48341-100ug

| | |
|---------------|--|
| G-102-C | Goat anti-Human IgG Secondary Antibody [Unconjugated] |
| NB7446 | Goat anti-Human IgG Fc Secondary Antibody |
| DDXCH01P-100 | Human IgG1 Isotype Control |
| NBP2-14363PEP | Bone marrow stromal cell antigen 1/CD157 Recombinant Protein Antigen |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-48341

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

