

# Product Datasheet

## XPR1 Antibody (HL2632) NBP3-48709

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-48709](http://www.novusbio.com/NBP3-48709)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-48709](http://www.novusbio.com/reviews/destination/NBP3-48709)



**NBP3-48709**

XPR1 Antibody (HL2632)


Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	HL2632
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS

Product Description	
Description	Novus Biologicals Rabbit XPR1 Antibody (HL2632) (NBP3-48709) is a recombinant monoclonal antibody validated for use in IHC, WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	9213
Gene Symbol	XPR1
Species	Human, Mouse, Rat
Immunogen	Recombinant protein encompassing a sequence within the Intracellular domain of human XPR1. The exact sequence is proprietary.


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:500-1:3000, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin

**Images**


Boiled and unboiled HCT116 whole cell (30 ug) were separated by 7.5% SDS-PAGE, and the membrane was blotted with XPR1 antibody [HL2632] (NBP3-48709) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.




Unboiled rat tissue extract (50 ug) was separated by 7.5% SDS-PAGE, and the membrane was blotted with XPR1 antibody [HL2632] (NBP3-48709) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody, and the signal was developed with Trident ECL plus-Enhanced.




Non-transfected (-) and transfected (+) unboiled 293T whole cell extracts (30 ug) were separated by 7.5% SDS-PAGE, and the membrane was blotted with XPR1 antibody [HL2632] (NBP3-48709) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.




XPR1 antibody [HL2632] detects XPR1 protein at cell membrane by immunohistochemical analysis. Sample: Paraffin-embedded mouse kidney. XPR1 stained by XPR1 antibody [HL2632] (NBP3-48709) diluted at 1:100. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



XPR1 antibody [HL2632] detects XPR1 protein at cell membrane by immunofluorescent analysis. Sample: HCT116 cells were fixed in 4% paraformaldehyde at RT for 15 min. Green: XPR1 stained by XPR1 antibody [HL2632] (NBP3-48709) diluted at 1:500. Blue: Fluoroshield with DAPI .




Boiled and unboiled mouse tissue extract (50 ug) were separated by 7.5% SDS-PAGE, and the membrane was blotted with XPR1 antibody [HL2632] (NBP3-48709) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody, and the signal was developed with Trident ECL plus-Enhanced.




Non-transfected (-) and transfected (+) 293T whole cell extracts (30 ug) were separated by 7.5% SDS-PAGE, and the membranes were blotted with XPR1 antibody [HL2632] (NBP3-48709) diluted at 1:1000 and competitor's antibody (Highly Cited Antibody) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.

\*The competitor is not affiliated with and does not endorse this product.



Boiled and unboiled mouse tissue extract (50 ug) were separated by 7.5% SDS-PAGE, and the membranes were blotted with XPR1 antibody [HL2632] (NBP3-48709) diluted at 1:1000 and the highly cited competitor's diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.

\*The competitor is not affiliated with and does not endorse this product.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-48709**

---

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-48709](http://www.novusbio.com/reviews/submit/NBP3-48709)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

