

Product Datasheet

P2Y2 Antibody (HL2801)

NBP3-48734

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-48734

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-48734



NBP3-48734

P2Y2 Antibody (HL2801)


Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	HL2801
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS

Product Description	
Description	Novus Biologicals Rabbit P2Y2 Antibody (HL2801) (NBP3-48734) is a recombinant monoclonal antibody validated for use in IHC and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	5029
Gene Symbol	P2RY2
Species	Human, Mouse
Immunogen	Synthetic peptide encompassing a sequence within the Extracellular domain mouse of P2Y2. The exact sequence is proprietary.


Product Application Details	
Applications	Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin 1:100-1:1000

Images


P2Y2 antibody [HL2801] detects P2Y2 protein by immunohistochemical analysis. Sample: Paraffin-embedded mouse tissues. P2Y2 stained by P2Y2 antibody [HL2801] (NBP3-48734) diluted at 1:200. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min




P2Y2 antibody [HL2801] detects P2Y2 protein by immunofluorescent analysis. Sample: Mock and transfected 293T cells were fixed in ice-cold MeOH for 5 min. Green: P2Y2 stained by P2Y2 antibody [HL2801] (NBP3-48734) diluted at 1:500. Blue: Fluoroshield with DAPI .




P2Y2 antibody [HL2801] detects P2Y2 protein by immunohistochemical analysis. Sample: Paraffin-embedded rat small intestine. P2Y2 stained by P2Y2 antibody [HL2801] (NBP3-48734) diluted at 1:200. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min




P2Y2 antibody [HL2801] detects P2Y2 protein by immunohistochemical analysis. Sample: Paraffin-embedded rat kidney. P2Y2 stained by P2Y2 antibody [HL2801] (NBP3-48734) diluted at 1:200. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min




P2Y2 antibody [HL2801] detects P2Y2 protein at cytoplasm by immunohistochemical analysis. Sample: Paraffin-embedded mouse eye. Green: P2Y2 stained by P2Y2 antibody [HL2801] (NBP3-48734) diluted at 1:200. Red: beta Tubulin 3/ Tuj1, a neural marker, stained by beta Tubulin 3/ Tuj1 antibody [GT11710] diluted at 1:500. Blue: Fluoroshield with DAPI .Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



P2Y2 antibody [HL2801] detects P2Y2 protein by immunofluorescent analysis. Sample: THP-1 cells were fixed in 4% paraformaldehyde at RT for 15 min. Green: P2Y2 stained by P2Y2 antibody [HL2801] (NBP3-48734) diluted at 1:100. Blue: Fluoroshield with DAPI .



P2Y2 antibody [HL2801] detects P2Y2 protein at cytoplasm by immunohistochemical analysis. Sample: Paraffin-embedded mouse hippocampus. Green: P2Y2 stained by P2Y2 antibody [HL2801] (NBP3-48734) diluted at 1:200. Red: beta Tubulin 3/ Tuj1, a neural marker, stained by beta Tubulin 3/ Tuj1 antibody [GT11710] diluted at 1:500. Blue: Fluoroshield with DAPI .Antigen Retrieval: Citrate buffer, pH 6.0, 15 min





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-48734

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
DYC1825-2	ERK1 [p Thr202, p Tyr204] [Biotin]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-48734

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

