

Product Datasheet

alpha-1D Adrenergic R/ADRA1D Antibody (HL3267) NBP3-48784

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-48784

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-48784



NBP3-48784

alpha-1D Adrenergic R/ADRA1D Antibody (HL3267)

Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	HL3267
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS


Product Description	
Description	Novus Biologicals Rabbit alpha-1D Adrenergic R/ADRA1D Antibody (HL3267) (NBP3-48784) is a recombinant monoclonal antibody validated for use in ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	146
Gene Symbol	ADRA1D
Species	Human
Immunogen	Synthetic peptide of human alpha-1D Adrenergic R/ADRA1D.

Product Application Details	
Applications	Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Immunocytochemistry/ Immunofluorescence 1:100-1:1000




Images

alpha-1D Adrenergic R/ADRA1D antibody [HL3267] detects alpha-1D Adrenergic R/ADRA1D protein by immunofluorescent analysis. Sample: A549 cells were fixed in ice-cold MeOH for 5 min. Green: alpha-1D Adrenergic R/ADRA1D stained by alpha-1D Adrenergic R/ADRA1D antibody [HL3267] (NBP3-48784) diluted at 1:500. Blue: Fluoroshield with DAPI .



alpha-1D Adrenergic R/ADRA1D antibody [HL3267] detects alpha-1D Adrenergic R/ADRA1D protein by immunofluorescent analysis. Sample: Mock and transfected 293T cells were fixed in ice-cold MeOH for 5 min. Green: alpha-1D Adrenergic R/ADRA1D stained by alpha-1D Adrenergic R/ADRA1D antibody [HL3267] (NBP3-48784) diluted at 1:500. Blue: Fluoroshield with DAPI .





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-48784

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-90229PEP	alpha-1D Adrenergic R/ADRA1D Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-48784

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

