

# Product Datasheet

## TNF RI/TNFRSF1A Antibody Pair [Biotin] NBP3-57712-1Set

Unit Size: 1 Set

Storage of components varies. See protocol for specific instructions.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-57712](http://www.novusbio.com/NBP3-57712)

Updated 9/29/2026 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-57712](http://www.novusbio.com/reviews/destination/NBP3-57712)



**NBP3-57712-1Set**

TNF RI/TNFRSF1A Antibody Pair [Biotin]

| Product Information  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
|----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size            | 1 Set                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Concentration        | Concentration of individual antibodies may be found on the vial label. If unlisted please contact technical services.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Storage              | Storage of components varies. See protocol for specific instructions.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Conjugate            | Biotin                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Product Description  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Description          | The kit contains a capture antibody, detection antibody, and recombinant protein standard for assay development. These antibodies have been validated to function together in an ELISA format using the provided standard. Optimal buffers, antibody concentrations, incubation times, temperatures, and assay conditions have not been predetermined, allowing investigators to optimize assay sensitivity and dynamic range based on their specific sample type and experimental objectives. The quantities of kit components are not matched, and individual components may be purchased separately. This product is for research use only. |
| Gene ID              | 7132                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Gene Symbol          | TNFRSF1A                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Species              | Guinea Pig                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Assay Type           | Sandwich ELISA                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Suitable Sample Type | Cell culture supernatants. With minor modifications to the assay procedure, this ELISA can also be used in additional biological sample types reported in the literature, including serum, plasma, and tissue homogenates, providing flexibility for a broad range of immunology, inflammation, and translational research applications.                                                                                                                                                                                                                                                                                                       |
| Sample Volume        | 100 ul                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Antibody Pairs are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-57712](http://www.novusbio.com/reviews/submit/NBP3-57712)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

