

# Product Datasheet

## CNPase Antibody (CL2887) - BSA Free NBP3-57807-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-57807](http://www.novusbio.com/NBP3-57807)

Updated 10/7/2026 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-57807](http://www.novusbio.com/reviews/destination/NBP3-57807)



**NBP3-57807-100ul**

CNPase Antibody (CL2887) - BSA Free

| Product Information         |                                                                                                                                       |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 100 ul                                                                                                                                |
| Concentration               | 1 mg/ml                                                                                                                               |
| Storage                     | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.                                                |
| Clonality                   | Monoclonal                                                                                                                            |
| Clone                       | CL2887                                                                                                                                |
| Preservative                | 0.02% Sodium Azide                                                                                                                    |
| Isotype                     | IgG2b                                                                                                                                 |
| Purity                      | Protein A purified                                                                                                                    |
| Buffer                      | PBS (pH 7.2), 40% glycerol                                                                                                            |
| Product Description         |                                                                                                                                       |
| Host                        | Mouse                                                                                                                                 |
| Gene ID                     | 1267                                                                                                                                  |
| Gene Symbol                 | CNP                                                                                                                                   |
| Species                     | Human, Mouse, Rat                                                                                                                     |
| Immunogen                   | This antibody was developed using a recombinant protein derived from P09543, with the exact immunogen sequence remaining proprietary. |
| Product Application Details |                                                                                                                                       |
| Applications                | Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry                                                                     |
| Recommended Dilutions       | Western Blot 1 ug/ml, Immunohistochemistry 1:10000 - 1:20000, Immunohistochemistry-Paraffin 1:10000 - 1:20000                         |
| Application Notes           | For IHC-Paraffin, HIER pH 6 retrieval is recommended.                                                                                 |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-57807](http://www.novusbio.com/reviews/submit/NBP3-57807)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

